
Article Number : C0311

■ **INTENDED USE** : The iF Helicobacter Test is intended to detect the Helicobacter pylori urease activity in biopsy specimens.

■ SUMMARY & PRINCIPLE OF THE TEST

The iF Helicobacter Test is an in vitro diagnostic test for detecting Helicobacter pylori in gastric mucosal biopsies. Helicobacter pylori is a spiral shaped Gram-negative bacterium that colonizes the human stomach. Infection with H. pylori causes chronic inflammation and significantly increases the risk of developing duodenal and gastric ulcer disease and gastric cancer. Infection with H. pylori is the strongest known risk factor for gastric cancer. H. pylori produces large amount of urease. If H. pylori is present in the biopsy specimens, urea would be converted to ammonia and carbon dioxide by the action of urease. The ammonia would cause pH changes which would be detected by a pH indicator present in the gel.

The iF Helicobacter Test is a cartridge with a single well prefilled with a selective agar gel containing urea, pH indicator, buffer and a bacteriostatic reagent preventing the growth of contaminating organisms.

■ MATERIALS SUPPLIED

The iF Helicobacter Test kit consists of 50 individual test cartridges and a package insert.

■ STORAGE

Store at 2~8°C ,DO NOT FREEZE!

The iF Helicobacter test is stable at 2~8°C until the expiration date (indicated on the package label) and at 25°C for 2 months.

■ SPECIMEN COLLECTION & PREPARATION

1. Wear disposable gloves when handling specimens and wash hands thoroughly afterwards. Collect specimens with biopsy forceps.
2. The usual area to biopsy is the sump of the antrum, along the greater curve.
3. The specimen size would be approximately 1-3mm. If the specimen appears to be very small, it may be worthwhile taking a second specimen.
4. Don't use the test after the expiration date.

■ PROCEDURE


1. Check if the gel is present in the well and its color is yellow.
2. The iF Helicobacter Test should be allowed to equilibrate to room temperature before using.
3. Peel back the label sticker from the plastic slide until you could see the yellow gel. DO NOT remove the label.

4. With a sterile needle, take the biopsy specimen from the biopsy forceps and push it into the gel.
5. Reseal the label sticker on the plastic slide by pressing so that the gel is covered.
6. Record the patient name, date and time on the label.
7. Incubate the plastic slide at 37°C for 2 hours for rapid results. Alternatively, incubate the test at room temperature for reading during 24 hours.
8. DO NOT read the results after 24 hours.

■ INTERPRETATION OF RESULTS


1. Positive Results

If H. pylori is present in the biopsy specimen, the yellow gel would be changed to the red or violet gel after incubation. The color change from yellow to red/violet indicates a positive result.

2. Negative Results

If H. pylori is absent in the biopsy specimen, the gel would remain yellow after incubation and keep yellow up to 24 hours at room temperature. This constitutes a negative result.

■ LIMITATION OF THE TEST

False Negative

The following conditions may lead to false negative results:

1. If patient has recently taken antibiotics, bismuth salts, proton pump inhibitors, or sucralfate which could inhibit the organism.
2. Very low numbers of H. pylori in the biopsy specimen.
3. Patchy H. pylori distribution so that the organism is not captured in the biopsy specimen.

False positive

The following conditions may lead to false positive results:

1. Patients with achlorhydria or pernicious anemia, previous gastric surgery or who recently has taken medicines (antacid, H2 receptor antagonist) may yield false positive results.
2. When acid is absent, commensal organisms such as Proteus spp., may grow in the stomach and produce urease. False positive results due to bacteria others than H. pylori would not react before two hours because they may produce urease much less than H. pylori.
3. If the specimen contains blood, it is usually seen as red edge around the specimen. This is not a positive result.

⚠ WARNING / NOTES

1. All reagents of this kit are strictly intended for in vitro diagnostic and professional use only.
2. The test is a qualitative test to detect H. pylori in biopsy specimens.
3. The test should be stored at 2~8°C.
4. Wear disposable gloves when handling specimens and wash hands thoroughly afterwards.
5. Do not eat, drink, smoke or apply makeup in areas where specimens or kit reagents are handled.
6. If the gel becomes pink BEFORE the specimen is added, the test should not be used.
7. DO NOT read the results after 24 hours when the specimen is added to the gel.
8. Disposal of medical wastes must be in accordance with all applicable federal, state, and local regulations.

■ REFERENCES

1. M. Castro Fernández Diagnosis of Helicobacter pylori infection using urease rapid test in patients with bleeding duodenal ulcer: influence of endoscopic signs and simultaneous corporal and antral biopsies. REV ESP ENFERM DIG 2004; 96(9): 599-605.
2. Lydia E. Wroblewski. Helicobacter pylori and Gastric Cancer: factors that modulate disease risk. CLINICAL MICROBIOLOGY REVIEWS. 2010;(4):713-739.

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