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1. M. Castro Fernández. Diagnosis of Helicobacter pylori infection using urease rapid test in patients with bleeding duodenal ulcer: influence of endoscopic signs and simultaneous corporal and antral biopsies. REV ESP ENFERM DIG 2004; 96(9): 599-605.
2. Lydia E. Wroblewski. Helicobacter pylori and Gastric Cancer: Factors That Modulate Disease Risk. CLINICAL MICROBIOLOGY REVIEWS. 2010;(4) :713-739.

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
Article number: C0811

Intended Use:

The IF Rapid Helicobacter Test is intended to detect the Helicobacter pylori urease activity in biopsy specimens.

1. Summary & Principle of The Test

The IF Rapid Helicobacter Test is an in vitro diagnostic test for detecting Helicobacter pylori in gastric mucosal biopsies. Helicobacter pylori is a spiral shaped Gram-negative bacterium that colonizes the human stomach. Infection with H. pylori causes chronic inflammation and significantly increases the risk of developing duodenal and gastric ulcer disease and gastric cancer. Infection with H. pylori is the strongest known risk factor for gastric cancer. H. pylori produces large amount of urease. If H. pylori is present in the biopsy specimens, urea would be converted to ammonia and carbon dioxide by the action of urease. The ammonia will cause pH changes which is detected by a pH indicator present in the reagent. The IF Rapid Helicobacter Test is a liquid reagent containing urea, pH indicator, buffer and a bacteriostatic reagent preventing the growth of contaminating organisms in the tube.

2. Materials Supplied

The IF Rapid Helicobacter Test kit includes of 50 individual test tubes and 50 fork-like sample collectors.

3. Storage

Store at 2 to 8°C ,DO NOT FREEZE!

The IF Rapid Helicobacter test is stable at 2 to 8°C until the expiration date (indicated on the package label).


4. Specimen Collection & Preparation

1. Wear disposable gloves when handling specimens and wash hands thoroughly afterwards. Collect specimens with biopsy forceps.
2. The usual area for biopsy is the sump of the antrum, along the greater curve of the stomach.
3. The specimen size would be approximately 1-3mm. If the specimen appears to be very small, it is advised to take a second specimen.
4. Do not use the test after the expiration date.

5. Procedure

1. Make sure that the reagent slightly exceeds the reaction area and that its color is yellow.
2. The IF Rapid Helicobacter Test should be allowed to equilibrate to room temperature before using.
3. Use sample collector to collect appropriate size of tissue and put it into the test tube.
4. Screw the cap tightly, the specially designed collector will center the specimen in the reagent and also push out any excess liquid.


5. Record the patient name, date and time on the sticker. Put the sticker on the tube while avoiding to cover the reagent zone.
6. Put the labeled sample tube on the rack which is provided in the reagent kit.
7. Interpret the result after 1 hour of incubation at room temperature. For easier reading check the results through the window of the rack.

6. Interpretation of Results


1. Positive Results

If H. pylori is present in the biopsy specimen, the yellow reagent would be changed to red or violet color after incubation. The color change from yellow to red/violet indicates a positive result.

2. Negative Results

If H. pylori is absent in the biopsy specimen, the reagent would remain yellow after incubation and keep yellow up to 24 hours at room temperature. This constitutes a negative result.

7. Limitation of The Test

False Negative

The following conditions may lead to false negative results:

1. If patient has recently taken antibiotics, bismuth salts, proton pump inhibitors, or sucralfate which could inhibit the organism.

2. Very low numbers of H. pylori in the biopsy specimen.
3. Patchy H. pylori distribution so that the organism is not captured in the biopsy specimen.

False positive

The following conditions may lead to false positive results:

1. Patients with achlorhydria or pernicious anemia, previous gastric surgery or who recently has taken medicines (antacid, H2 receptor antagonist) may yield false positive results.
2. When acid is absent, commensal organisms such as Proteus spp., may grow in the stomach and produce urease. False positive results due to bacteria others than H. pylori would not react before two hours because they may produce urease much less than H. pylori.
3. If the specimen contains blood, it is usually seen as red edge around the specimen. This is not a positive result.

⚠ Precautions / Notes

1. All reagents of this kit are strictly intended for in vitro diagnostic and professional use only.
2. The test is a qualitative test to detect H. pylori in biopsy specimens.
3. The test should be stored at 2 to 8°C.
4. Wear disposable gloves when handling specimens and wash hands thoroughly afterwards.
5. Do not eat, drink, smoke or apply makeup in areas where specimens or kit reagents are handled.
6. If the reagent becomes pink BEFORE the specimen is added, the test should not be used.
7. DO NOT read the results more than 24 hours after the specimen was added to the reagent.
8. Disposal of medical wastes must be in accordance with all applicable federal, state, and local regulations.
9. In the most cases, positive results will already be seen at 5 and 30 minutes. All positive results are shown at 1 hour.

References

1. M. Castro Fernández Diagnosis of Helicobacter pylori infection using urease rapid test in patients with bleeding duodenal ulcer: influence of endoscopic signs and simultaneous corporal and antral biopsies. REV ESP ENFERM DIG 2004; 96(9): 599-605.
2. Lydia E. Wroblewski. Helicobacter pylori and Gastric Cancer: factors that modulate disease risk. CLINICAL MICROBIOLOGY REVIEWS. 2010;(4):713-739.

■ Revision: 1st edition

IN FUNG Co., Ltd.

No. 61-1, 40th Rd., Taichung Industrial Park, Taichung City 40768, Taiwan R.O.C.

Tel: +886-(0)4-23550886

Fax: +886-(0)4-23550996

Web: www.infung.com.tw