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1. M. Castro Fernández. Diagnosis of Helicobacter pylori infection using urease rapid test in patients with bleeding duodenal ulcer: influence of endoscopic signs and simultaneous corporal and antral biopsies. REV ESP ENFERM DIG 2004; 96(9): 599-605.
2. Lydia E. Wroblewski. Helicobacter pylori and Gastric Cancer: Factors That Modulate Disease Risk. CLINICAL MICROBIOLOGY REVIEWS. 2010;(4):713-739.

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"IF" Rapid Helicobacter Disc Test (Non-Sterile)


Intended Use:

The "IF" Rapid Helicobacter Disc Test (RHD Test) is intended to detect the Helicobacter pylori urease activity in biopsy specimens.

1. Summary & Principle of The Test

The "IF" RHD Test is an in vitro diagnostic test for detecting Helicobacter pylori in gastric mucosal biopsies. Helicobacter pylori is a spiral shaped Gram-negative bacterium that colonizes the human stomach. Infection with H. pylori causes chronic inflammation and significantly increases the risk of developing duodenal and gastric ulcer disease, and gastric cancer. Infection with H. pylori is the strongest known risk factor associated with gastric cancer. H. pylori produces large amount of urease. If H. pylori is present in the biopsy specimens, urea would be converted to ammonia and carbon dioxide by the action of urease. The ammonia would cause pH changes which would be detected by a pH indicator present in the RHD Test.

2. Specification

Product number	C0914	C0911	C0926
specification	Rapid Helicobacter Disc Test		
	35 T	50T	70 T

3. Storage

Store at 2 to 30°C, DO NOT FREEZE!

The RHD Test is stable at 2 to 30°C until the expiration date (indicated on the package label).

4. Specimen Collection & Preparation

1. Wear disposable gloves when handling specimens and wash hands thoroughly afterwards. Collect specimens with biopsy forceps.
2. The usual area to biopsy is the sump of the antrum, along the greater curve.
3. The specimen size should be approximately 1-3mm. If the specimen appears to be very small, a second specimen is highly recommended.
4. Do not use the test after the expiration date.

5. Procedure

1. Check if the test disc is intact in the well and its color is yellow.
2. The RHD Test should be allowed to equilibrate to room temperature before use.
3. Peel back the label sticker from the plastic slide until you could see the yellow test disc in the well. Do NOT remove the plastic disc from the sticker.
4. Put the specimen into the well and make sure that the specimen and the test disc are fully contact.
5. Reseal the label sticker and press the plastic disc to squeeze out the tissue fluid of the specimen.
6. Record the patient name, date and time on the label.
7. Interpret the result after 1 hour of incubation at room temperature.
8. DO NOT read the results after 2 hours.

6. Interpretation of Results

1. Positive Results

If H. pylori is present in the biopsy specimen, the yellow color test disc would be changed to red or violet color after incubation. This color change from yellow to red/violet indicates a positive result.

2. Negative Results

If H. pylori is absent in the biopsy specimen, the test disc would remain yellow after incubation and keep yellow up to 2 hours at room temperature. This indicates a negative result.

7. Limitation of The Test

False Negative

The following conditions may lead to false negative results:

1. If patient has recently taken antibiotics, bismuth salts, proton pump inhibitors, or sucralfate, which could inhibit the organism.
2. Very low number of H. pylori organism in the biopsy specimen.
3. Patchy H. pylori distribution so that the organism is not captured in the biopsy specimen.

False positive

The following conditions may lead to false positive results:

1. Patients with achlorhydria or pernicious anemia, previous gastric surgery or who recently has taken medicines (antacid, H2 receptor antagonist) may yield false positive results.
2. When acid is absent, commensal organisms such as Proteus spp., may grow in the stomach and produce urease. False positive results due to bacteria other than H. pylori may produce urease. But the amount of urease produced by these bacteria are much lower than H. pylori and will not cause a yellow to red/violet color change within two hours.
3. If the specimen contains blood, it is usually seem as red edge around the specimen. This is not a positive result.

⚠ Precautions / Notes

1. All tests of this kit are strictly intended for in vitro diagnostic and professional use only.
2. The test is a qualitative test to detect H. pylori in biopsy specimens.
3. The test should be stored at 2 to 30°C.
4. Wear disposable gloves when handling specimens and wash hands thoroughly afterwards.
5. Do not eat, drink, smoke or apply makeup in areas where specimens or tests are handled.
6. If the test disc becomes pink BEFORE the specimen is added, the test should not be used.
7. DO NOT read the results more than 2 hours after the specimen was added to the test disc.
8. Disposal of medical wastes must be in accordance with all applicable federal, state, and local regulations.

References

1. M. Castro Fernández Diagnosis of Helicobacter pylori infection using urease rapid test in patients with bleeding duodenal ulcer: influence of endoscopic signs and simultaneous corporal and antral biopsies. REV ESP ENFERM DIG 2004; 96(9): 599-605.
2. Lydia E. Wroblewski. Helicobacter pylori and Gastric Cancer: factors that modulate disease risk. CLINICAL MICROBIOLOGY REVIEWS. 2010;(4):713-739.

■ Revision: Second edition (Oct 2024)

IN FUNG Co., Ltd.

No. 61-1, 40th Rd., Taichung Industrial Park, Taichung City 40768, Taiwan R.O.C.

Tel: +886-(0)4-23550886

Fax: +886-(0)4-23550996

Web: www.infung.com.tw